



# MAcURAY

## Specification

Power . . . . .	200~240V AC 50/60Hz
Current Consumption . . . . .	Cooling Device : 150W
. . . . .	Smoke Evacuator : 1,100W
Air Cooling Delivery . . . . .	45 CFM
Vacuum Power . . . . .	2,400 mmH2O
Range of Temperature . . . . .	Air Cooling : Max -14°C
. . . . .	Contact Probe : Cooling Max -29°C
. . . . .	: Warm +35 ~ +43°C
. . . . .	(based on +22°C Surrounding Temp.)
. . . . .	(Leading Setting Time: 10sec. ~ 30minutes)
Weight . . . . .	42 kg(Net.)
Dimension . . . . .	320(W) x 350(D) x 990(H) / mm
. . . . .	1290(H) mm With Pole
. . . . .	(Available Separation Pole)



### Multi Postcare System

- Air Cooling
- Cooling Contact Probe
- Smoke Evacuator



# MAcURAY

Electrical Cold Therapy Device

## Cooling Contact Probe

Quick Down & Low Temperature  $-29^{\circ}\text{C}$   
Wide rubbing SUS Contact Probe

## Air Cooling

Continuously Keep under  $-14^{\circ}\text{C}$   
Strong Air Outlet 45 CFM  
Various Air Outlet Modules  
No Noise

## Smoke Evacuator

Easy On/Off Control by Foot Switch  
Strong Vacuum & Minimized Noise  
Adopt a High Performance Smoke Filter  
Adoptable Wide Smoke Capture Module

# Cooling by TEC Module & Environment Friendly

- Simple & Free angle adjustment of the Air cooling outlet.
- Air cooling & contact type probe can use simultaneously.
- Always supply clean air by adopting inside air filter.
- Convenient monitoring and controlling with full touch LCD Screen.
- Short standby time to operate.
- Strong smoke evacuator by adopting a high performance smoke filter.
- Multipurpose device



Three in One System...

Simply

Easily

## What is MAcURAY?

MAcURAY is a Unique Electrical Cold Therapy Device including Smoke Evacuator.

It reduces Pain by maintaining low temperature with Air Cooling & Contact Typed Probe

MAcURAY Smoke Evacuator can inhale smoke occurs during Laser & IPL treatment. Furthermore, Air cooling & cooling contact probe can be used to calm down hot skin after skin treatment like Laser, IPL, derma roller and stamp.

As full touch LCD screen adopted, you can easily select functions and control temperature

## MAcURAY Multi Postcare System

### 3 functions Can be . . . . . used simultaneously

### MAcURAY Combine with Laser/IPL System

Minimized Edema, Redness, Thermal Injury, and Reduce Pain.  
Air Cooling & Cooling Contact Probe help Rapid Recovery from wound parts.

### MAcURAY Various Application

Applicable to any part of operation  
Relief Pain, Reduce Skin Stimulation & Edema

